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1. USP <645> Water Conductivity. Retrieved February 14, 2019 from <https://hmc.usp.org/sites/default/files/documents/HMC/GCs-Pdfs/c645.pdf>
2. Sievers Lean Lab: Simultaneous Stage 1 Conductivity and TOC Lab Testing of Pharmaceutical Water (300 40030). Retrieved February 14, 2019 from [https://www.suezwatertechnologies.com/kcpguest/documents/Application%20Notes\\_Cust/Americas/English/ANai\\_300\\_40030\\_EN.pdf](https://www.suezwatertechnologies.com/kcpguest/documents/Application%20Notes_Cust/Americas/English/ANai_300_40030_EN.pdf)
3. DUCT Vial Performance and Stability (300 00297). Retrieved February 14, 2019 from [https://www.suezwatertechnologies.com/kcpguest/documents/Technical%20Bulletins\\_Cust/Americas/English/TBai\\_300\\_00297\\_EN.pdf](https://www.suezwatertechnologies.com/kcpguest/documents/Technical%20Bulletins_Cust/Americas/English/TBai_300_00297_EN.pdf)
4. Reserve Sample Bottles for Conductivity and TOC (300 00299). Retrieved February 14, 2019 from [https://www.suezwatertechnologies.com/kcpguest/documents/Technical%20Bulletins\\_Cust/Americas/English/TBai\\_300\\_00299\\_EN.pdf](https://www.suezwatertechnologies.com/kcpguest/documents/Technical%20Bulletins_Cust/Americas/English/TBai_300_00299_EN.pdf)
5. Sievers Inorganic Carbon Remover (ICR) (300 00109). Retrieved February 14, 2019 from [https://www.suezwatertechnologies.com/kcpguest/documents/Application%20Notes\\_Cust/Americas/English/ANai\\_300\\_00109\\_EN.pdf](https://www.suezwatertechnologies.com/kcpguest/documents/Application%20Notes_Cust/Americas/English/ANai_300_00109_EN.pdf)
6. Low Level Linearity Conductivity Study on the Sievers M9 TOC Analyzer (300 00339). Retrieved February 14, 2019 from [https://www.suezwatertechnologies.com/kcpguest/documents/Application%20Notes\\_Cust/Americas/English/ANai\\_300\\_00339\\_EN.pdf](https://www.suezwatertechnologies.com/kcpguest/documents/Application%20Notes_Cust/Americas/English/ANai_300_00339_EN.pdf)
7. UPW 07-10 Multiple Products Biological Contamination (800 19025) Retrieved February 14, 2019. [https://www.suezwatertechnologies.com/kcpguest/documents/Technical%20Bulletins\\_Cust/Americas/English/ai\\_UPW\\_07-10\\_EN.pdf](https://www.suezwatertechnologies.com/kcpguest/documents/Technical%20Bulletins_Cust/Americas/English/ai_UPW_07-10_EN.pdf)

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